

Longest™

Compression Therapy System DVT Sleeve



Longest™ Compression Therapy System DVT Sleeve (single patient use) as an essential part of Longest™ Compression Therapy System is compatible with both Longest™ Compression Therapy Device LGT-2200DVT and LGT-2201DVT.

It is made of premium fabric and TPU with featuring lightweight, breathable, comfortable which is able to help prevent a build-up of heat and moisture contributing to increase patient's compliance during therapy.

With an extensive range of foot, calf, thigh, foot and calf, thigh and calf sleeves, Longest™ DVT Sleeve address a wide range of patient types, and caregiver needs and preferences.

⊕ Feature

- ⊗ Made of premium fabric with more lightweight, breathable and comfortable
- ⊗ Compatible with both Longest™ LGT-2200DVT and LGT-2201DVT
- ⊗ With a wide range of types for different application
- ⊗ One-Three (1-3) chambered sleeve with intermittent and sequential (gradient) inflation

DVT Sleeve Size Details

Order No.	Sleeve	Size	Proximal Circumference
220092	Thigh Length Sleeve	S	<55.9cm
220092M		M	<71.1cm
220092L		L	<91.4cm
220093M	Calf and Foot Sleeve	M	<53.3cm
220094	Calf Sleeve	S	<40.6cm
220094M		M	<48.3cm
220094L		L	<55.9cm
220095M	Foot Sleeve	M	<40.6cm



Longest
PHYSIO & REHAB DEVICES

Guangzhou Longest Medical Technology Co., Ltd.

Address: 301 & 401 of Building 2 & Building 3, No.96, Chuangqiang Road, Ningxi Street, Zengcheng District, Guangzhou, Guangdong Province, 511399, China

Tel: 86-020-66353999

Fax: 86-020-66353920

E-mail: export@longest.cn

Website: www.longestmedical.com

